

Ladd Research  
3 Ewing Place  
Essex Junction, VT 05452  
Tel: (802) 658-4961  
Email: sales@laddresearch.com  
Web: www.laddresearch.com



## Catalog Number 55140: Durcupan ACM Embedding Kit Instructions for Use

These instructions were extracted from information supplied by Sigma/Fluka, St. Louis, MO

Mixtures of Durcupan ACM present an interesting alternative to methacrylates for embedding sections for electron microscopy. Durcupan hardens uniformly, practically without shrinkage. All fixing agents normally used for electron microscopic work are also suitable for the Durcupan ACM embedding process.

The embedding process consists of 3 phases: dehydration of the tissue with acetone or ethanol, infiltration with the embedding agent and curing.

The following 2 solutions should be prepared immediately prior to embedding:

### Solution 1

Durcupan ACM Resin - 10ml  
DDSA - 10ml  
Dibutyl phthalate - 0.1 to 0.2ml

### Solution 2

Durcupan ACM Resin - 10ml  
DDSA - 10ml  
Dibutyl phthalate - 0.1 to 0.2ml  
DMP 30 - 0.3 to 0.4ml

### Dehydration

1. 30% acetone – 15 minutes
2. 50% acetone – 30 minutes
3. 70% acetone – 30 minutes
4. 90% acetone – 30 minutes
5. Anhydrous acetone – 30 minutes. (Ethanol may be used)
6. Anhydrous acetone – 30 minutes

### Infiltration

7. **Three** parts anhydrous acetone (or ethanol) in **one** part Durcupan Solution 1. Room temperature. One hour.
8. **Two** parts anhydrous acetone (or ethanol) in **two** parts Durcupan Solution 1. Room temperature. One hour.
9. **One** part anhydrous acetone (or ethanol) in **three** parts Durcupan Solution 1. Room temperature. One hour. This period may be extended overnight if you wish.
10. Durcupan Solution 1. 50° C. 1-2 hours.

11. Durcupan Solution 1. 50° C. 1-2 hours.
12. Durcupan Solution 2. 50° C. 1-2 hours.

### **Embedding**

13. Remove the tissue from step 12 and place it in a dry gelatin or BEEM capsule.
14. The capsule is then filled with Durcupan Solution 2. After closing the capsule, harden it in a drying oven at 50-80° C for at least 48 hours.

In this manner blocks of a pale golden color are obtained with a hardness similar to that of the methacrylates. The cutting properties are also superb.

### **NOTES**

1. Durcupan Solutions 1 and 2 are preferably prepared at least 15 minutes before use and kept in a drying oven at 50° C.
2. Glassware used for measuring and mixing should be placed in acetone or absolute ethanol immediately after use and cleaned.
3. Durcupan mixtures are more soluble in acetone than in ethanol.
4. Addition of dibutyl phthalate reduces the brittleness of blocks and improves their cutting properties.
5. With materials containing minimal water, it is permissible to begin with step 2 (on the previous page).

### **CAUTION!!**

Take great care when working with Durcupan ACM. Do not breathe the vapor and avoid skin contact, because this may cause skin irritation and allergic reactions. Splashes on the skin must be washed immediately with a 3% boric acid solution. Frequent washing of hands, arms and face with lukewarm soapy water is recommended.