

Ladd Research
3 Ewing Place
Essex Junction, VT 05452
Tel: (802) 658-4961
Email: sales@laddresearch.com
Web: www.laddresearch.com



SEM Performance Standard Instructions

Catalog Number: 53130

1. Magnification calibration is straightforward. Specimen dimensions are specified in microns so magnification can be measured directly from photomicrographs in millimeters. Keep in mind the possibility of camera reduction ratios. It will be necessary to use large spot size for good definition at low magnifications.
2. The layout of this device is easy to remember. Start by using the sales brochure as a map. Simple stage maneuvers will allow the next range of interest to be placed in the field of view quickly.
3. The eight-point star in the upper right-hand corner of the specimen makes it possible to visualize where magnification center and scan rotation center occurs.
4. Pincushion or barrel distortion can easily be visualized by selecting any of the grid areas. The squares in these areas have reasonable edges. A straight edge placed upon a micrograph should easily touch all grid boundaries within the accuracy of the etching process. Any evidence of bowing would be considered distortion. Checking at various magnifications will verify if these conditions exist in the optics. Measuring along the X and Y or the X plus Y-axis will show geometry errors.
5. Depth of field measurements can be done with any of the grids; however, the ten-micron grid may be most useful. Tilting the specimen to 89° and looking across the bars will give calibrated measurements away from focus.
6. Tilt correction is best done using the circles, since perfect correction is not available, circles will allow you to visualize the best position of scan rotation vs. tilt correction for give KV and working distance.
7. Dynamic focus performance is best done in conjunction with tilt correction calibration. This will assure Y scan alignment along the vertical axis.
8. Orthogonality and linearity tests use the eight-point star. The intersecting angles of this device are 45° . They should remain that regardless of scan rotation or scan speed settings. Any deviation between the 90° intersection of the X and Y lines on the star should be interpreted as orthogonality errors usually due to the column scan coils, the display scan coil or an electronic problem in the deflection amplifiers or their power supplies.
9. Stage motion tests are accomplished using either the five-point star or for more critical tests use an appropriate set of cross hairs.
10. Video performance can be measured by placing the 5-micron squares in the center of the field of view and reducing the magnification until all detail on this target disappears. This will be the lower limit of

